

FP 431

Features:

- Closely controlled particle size
- Auto-polymerises with the recommended liquid
- Available in many different shades

Benefits:

- Cold Cure acrylic denture resins are believed to provide an effective means of repairing fractured dentures
- Highly accurate colour matching

Typical Properties

Makevale® Group Grades	FP 431
Median Particle Size (μ)	42 - 51
Residual Peroxide %	1.75 - 2.5
Solution Viscosity (Poise)	0.25 - 0.35
Recommended Liquid	Cold Cure Liquid

