

## FP 541

### Features:

- Closely controlled particle size
- Can be processed using either long or short curing techniques in boiling water

### Benefits:

- Heat Cure acrylic denture resins are designed for the manufacture of full or partial dentures
- Highly accurate colour matching

### Typical Properties

Makevale® Group Grades	FP 541
Median Particle Size ( $\mu$ )	40 - 60
Residual Peroxide %	0.25 - 0.38
Solution Viscosity (Poise)	2 - 3
Recommended Liquid	Heat Cure Denture base or Heat Cure Clear Liquid

